

Sampling and viral RNA extraction kit

TAAG X11 IVD

- **TAAG X11 IVD** is available in 100 tubes format, and includes sampling swabs

CAT: TAGC01007 - 100 Reactions

- **Technical Specifications**

| | |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product composition | 2 mL of transport buffer plus TAAG-Q Beads |
| Intended use | Stabilization, transport, deactivation and RNA extraction of clinical samples for SARS-CoV-2 detection Fluorescent probes |
| Type of clinical sample | Upper respiratory specimens such as Nasopharyngeal, Oropharyngeal, Nasal swabs and Saliva samples |

- Keep between 2 - 25 °C until expiration date

- **Research Use Only**

- **United States**

3710 Illinois Avenue, Unit A, St. Charles, IL.
60174 United States

Phone: +1 630 506.61.29

e-mail: contact_US@taag-genetics.com

